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## Enterobacteriaceae Count Plate

Type of Product	Product No.	Packaging Specifications
Enterobacteriaceae Count Plate®	HP013	20 plates/pack

### 【Usage】

Handy Plate® Enterobacteriaceae count Plates are pre-prepared, ready-to-use culture media products containing standard culture medium, cold-water-soluble gelling agent, pH indicator and colorimetric indicator. This product can be used for enumeration of Enterobacteriaceae in food and beverages.

### 【Instruction】

#### 1. Sample Preparation

Place 25 g (mL) of sample into a sterile homogenizing cup or bag containing 225 mL of buffered peptone water (or phosphate buffer, physiological saline). Homogenize for 1-2 min to prepare a 1:10 sample solution. Adjust the pH of the sample solution to 6.8-7.2 with sterile 1 mol/L NaOH or 1 mol/L HCl solution If needed. Use a 1 mL sterile pipette or tips to transfer 1 mL of the 1:10 sample solution to a test tube containing 9 mL of diluent. Shake to obtain a 1:100 sample solution. Prepare 10-fold serial dilutions of the sample solution. Replace the pipette after each dilution.

#### 2. Inoculation

Based on the estimated contamination level of the sample, select 2-3 dilutions for testing. Place a *Enterobacteriaceae* count plate horizontally, peel off the upper film, and use a sterile pipette to drop 1 mL of the sample solution onto the center of the plate. Gently lower the upper film to allow the sample solution to spread and reform into a gel and press lightly to expel bubbles.

(A few bubbles at the edges will not affect the results. Wait until the liquid fully diffuses before moving the plate).

#### 3. Incubation


Incubate plates horizontally with front side up at  $36 \pm 1^\circ\text{C}$  for 18 – 24 h. Do not stack more than 20 plates.

#### 4. Interpretation


Count red colonies with gas bubbles and yellow acid rings, regardless of the size and clarity. The optimal counting range is 15–150 CFU.

If colony count >150 CFU, estimate by counting representative grids and multiply by 20. If the whole plate is yellow with tiny red dots or bubbles, further dilute the sample solution. For further analysis of isolated colonies, please pick up the colonies for use.

### 【Process Diagram】


1. Cut along the dotted line with scissors and take out Handy plate.


2. Place a Handy plate horizontally, peel off the upper film, and drop 1 mL of sample solution to the plate.

3. Gently lower the upper film and press lightly to expel bubbles.

4. Incubate plates horizontally with front side up at  $36 \pm 1^\circ\text{C}$  for 18 – 24 h. Do not stack more than 20 plates.

5. Determine whether the target organism has grown on Handy plate according to the interpretation manual, and count the colonies accordingly.

Note: Scan the QR code on the left to watch the operation video.

### 【Storage and Shelf Life】

Store sealed at 2-8°C, valid for 18 months.

Put the unused Handy plates back into the packaging bag, fold it and seal it with a self-adhesive sticker. Once opened, use as soon as possible. Storage time shall not exceed 4 weeks.

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Handy plate may contain microorganisms after use and must be autoclaved by high pressure steam at 121°C for 30 min and then handled in strict accordance with regulations.